


THE T-BIO CONTROLLER

Automation controller box designed to regulate the temperature via a silicone band for the CellReady Bioreactor Mobius 3L. Data logging with a double PT100 probe and connection to Applikon or Sartorius process control unit.





TEMPERATURE CONTROLLER

Power supply	90 à 260 V AC
PID	Yes, Auto Adaptive
Protection	IP 50
Output	PT 100 for data logging
Sensor	
Temperature	2 x PT100
Heating Element	
Structure	Silicone ; 150W
Dimmensions	
Length (mm)	110
Depth (mm)	200
Height (mm)	170